A SWISS BRAND

# AUTOMATIC BLOOD PRESSURE MONITOR

MODEL M-400A  
USER'S MANUAL

[www.dr-frei.com](http://www.dr-frei.com)

**EN AUTOMATIC BLOOD PRESSURE MONITOR**

Model M-400A

User's manual (4-27)

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## DEAR DR. FREI AUTOMATIC BLOOD PRESSURE MONITOR OWNER

Thank you for choosing Automatic Blood Pressure Monitor **Dr. Frei**® TM model M-400A. We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Swiss Trademark **Dr. Frei**®.

This system is a fully automatic, digital blood-pressure measuring device for use on the upper arm, it enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measuring. This monitor features Blood Pressure Level Indicator, CAM technology - Clinically Accurate Measurement, Irregular Heartbeat Detection Advanced, alarm system (2 reminder sounds) and function of average of day/night/all stored results as well as 90 memory cells.

Before starting to use this device please study the user's manual carefully. The user's manual offers all information you need to measure your blood pressure and pulse correctly. For all questions concerning the device please contact your local distributor or **Dr. Frei**® service center in your country.

### ATTENTION

This blood pressure monitor is designed to carry out self-control over blood pressure but NOT to make self-diagnosis of hypertension/hypotension. Please DO NOT diagnose by yourselves basing on the measurement results obtained with the blood pressure monitor. Please DO NOT execute self-treatment of high/low blood pressure and DO NOT change the methods prescribed without consulting your doctor.

Type BF applied part.

Read the instructions carefully before using this device.

# TABLE OF CONTENTS

---

<b>IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT .....</b>	<b>6</b>
What Is Blood Pressure? .....	6
Which Values Are Normal? .....	7
<b>ADVANTAGES OF AUTOMATIC BLOOD PRESSURE MONITOR MODEL M-400A .....</b>	<b>9</b>
Blood Pressure Level Indicator .....	9
CAM technology - Clinically Accurate Measurement .....	9
Irregular Heartbeat Detection Advanced .....	9
<b>GETTING READY FOR MEASUREMENT .....</b>	<b>10</b>
Safety Precautions .....	10
Description of the Blood Pressure Monitor .....	11
Inserting the Batteries .....	12
Using a Mains Adaptor .....	12
Setting the Date and Time .....	13
Setting the Alarm .....	14
<b>MEASUREMENT PROCEDURE .....</b>	<b>15</b>
Fitting the Cuff .....	15
Taking a Single Measurement .....	16
CAM (Clinically Accurate Measurement) Mode .....	18
<b>MEMORY FUNCTION .....</b>	<b>20</b>
Average of All Stored Memories .....	20
Average of All Stored Day-Time Memories .....	20
Average of All Stored Night-Time Memories .....	21
Viewing the Stored Values .....	21
Delete Memories .....	22
<b>ERROR MESSAGES/TROUBLESHOOTING .....</b>	<b>22</b>
<b>CARE AND MAINTENANCE .....</b>	<b>24</b>
<b>SYMBOL INFORMATION .....</b>	<b>25</b>
<b>TECHNICAL SPECIFICATIONS .....</b>	<b>26</b>
<b>WARRANTY .....</b>	<b>27</b>


# IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT

## What Is Blood Pressure?


Blood pressure is the pressure that your blood extends to the vascular walls. Blood pressure is necessary to provide for constant blood flow inside the body. Thanks to it the cells get oxygen that provides for their normal functioning. The heart performs the function of a «pump», sending blood to the blood vessels. Each heart beat creates a certain level of the blood pressure.

There are 2 kinds of blood pressure: a systolic (upper) one, which corresponds to the heartbeat pushing blood into the arteriae; and a diastolic (lower) one, which means the blood pressure between two heartbeats.

Blood pressure is subject to fluctuations during the day even in healthy people. The fluctuations are influenced by a number of factors - time of day, person's condition, physical or mental activity, environment, etc.

## Day-Night Fluctuations of the Blood Pressure

An increase of blood pressure increases the burden onto the heart, affects blood vessels making their walls thick and less elastic.

One of the features of the hypertension is its ability to remain unnoticed for the patient at its early stages. That's why the self control of the blood pressure is so important. With the illness progressing, headaches and regular dizziness appear, the sight declines, the functioning of vitalses (encephalon, heart, kidneys, blood vessels) breaks down. Without special treatment the complications of hypertension might be kidney damages, breast-pang, paralytic stroke, aphasia, dementia, heart attack and stroke.

## Which Values are Normal?

Monitoring your routine blood pressure trend helps you to know your body condition. Human blood pressure naturally increases after reaching middle age.

This symptom is a result of continuous ageing of the blood vessels. Further causes include obesity, lack of exercise and cholesterol (LDL) adhering to the blood vessels, diabetes. Rising blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction.

Definitions and Classification of blood pressure levels according to 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension:. Data in mmHg.

Category	Systolic (mmHg)		Diastolic (mmHg)
Optimal	< 120	and	< 80
Normal*3	120-129	and/or	80-84
High normal	130-139	and/or	85-89
Grade 1 hypertension	140-159	and/or	90-99
Grade 2 hypertension	160-179	and/or	100-109
Grade 3 hypertension	≥ 180	and/or	≥ 110

Isolated systolic hypertension	$\geq 140$	and	$< 90$
--------------------------------	------------	-----	--------

Isolated systolic hypertension should be graded (1, 2, 3) according to systolic blood pressure values in the ranges indicated, provided that diastolic values are  $< 90$ mmHg.

Source: The European Society of Hypertension and European Society of Cardiology Task Force Members. 2007 ESH-ESC Practice Guidelines for the Management of Arterial Hypertension. *J Hypertens* 2007; 25: 1751-1762.

- The diagnosis of hypertension requires from the patient to combine medical treatment prescribed by the doctor and mode of life correction.
- People with normal pressure and high normal pressure are recommended to carry out self-control of their tension in order to timely take measures to decrease the blood pressure level down to the optimal one without using any medications.
- For people more than 50 years old high level of systolic blood pressure (higher than 140 mmHg) is more crucial than diastolic pressure.
- Even with blood pressure being normal, people run the bigger risk of hypertension development with advancing age.

#### **ATTENTION**

If you have normal results of blood pressure measured under calm conditions but your results are excessively high when measured under the conditions of physical or mental exhaustion, this might be a sign of so called brittle (that is unstable) hypertension. If you suspect that, please consult your doctor.

When measured correctly, if diastolic blood pressure is more than 120 mmHg, it is necessary to call the doctor immediately.

# ADVANTAGES OF AUTOMATIC BLOOD PRESSURE MONITOR MODEL M-400A

---

## Blood Pressure Level Indicator

---

Blood pressure level indicator is located along the left side in the display. The classification corresponds to ESH (European Society of Hypertension) ranges described in the table of the section «Which values are normal?». After the measurement there appears the color bars in the left part of the display: green zone - optimal blood pressure, yellow - elevated, orange - too high, red - dangerously high. This function helps you to self-orient in the measurement results.

## CAM technology - Clinically Accurate Measurement

---

CAM Technology - Clinically Accurate Measurement - provides you with the maximally accurate results of blood pressure measurement at home. Cardiologists recommend always take several repeated measurements and calculate their average values. The device in the CAM mode automatically takes 2 measurements in succession, analyzes the results and displays the ready values.

## Irregular Heartbeat Detection Advanced

---

This function indicates allorhythmic heartbeating. If the symbol of IHB Advanced  appears on the display that means that certain abnormality in heart beat frequency was detected during the measurement. In this case, the result may deviate from your normal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol  appears on a regular basis (e.g. several times a week with measurements taken daily) we advise you to tell your doctor. Please show your doctor the following explanation:

**Information for the doctor on frequent appearance of the Irregular**

### **Heartbeat Detector.**

This instrument is an oscillometric blood pressure monitor that also analyses pulse frequency during measurement. The instrument is clinically tested. The symbol  is displayed after the measurement, if pulse irregularities occur during measurement. If the symbol  appears more frequently (e.g. several times per week on measurements performed daily) we recommend the patient to seek medical advice. The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

## **GETTING READY FOR MEASUREMENT**

---


### **Safety Precautions**

---

1. Use this device ONLY for the intended use described in this manual.
2. Do NOT use accessories which are not specified by the manufacturer.
3. Do NOT use the device if it is not working properly or damaged.
4. Do NOT use the equipment where aerosol sprays are being used, or where oxygen is being administered.
5. Do NOT use under any circumstances on newborns.
6. This device does NOT serve as a cure for any symptoms or diseases. The data measured are for reference only. Always consult your doctor to have the results interpreted.
7. Keep the equipment and its flexible cord away from hot surfaces.
8. Do NOT apply the cuff to areas other than the place directed.
9. Proper maintenance and periodically calibration are essential to the longevity of your device. If you are concerned about your accuracy of measurement, please contact local customer service for help.
10. This system is intended for individuals age 18 or above.
11. Do not use it for diagnosis of hypertension or testing on newborns, babies, young children or persons who cannot express their consent.
12. If you have been diagnosed with a severe arrhythmia or irregular heartbeat, atrial or ventricular premature beats or atrial fibrillation, measurements made with this instrument should only be evaluated after consultation with the doctor or your healthcare professionals.

KEEP THIS MANUAL DURING THE WHOLE LIFETIME OF THE DEVICE


## Description of the Blood Pressure Monitor

- ① Cuff socket
- ② M/SET button
- ③ CAM button
- ④ START button
- ⑤ Battery cover
- ⑥ DC jack
- ⑦ LCD display


### Cuff:

M-400A offers a wide range soft conical cuff with D-ring, measuring size: 22~38cm, suitable for most of arm circumferences.

### Display screen:

- ⑧ Time Date
- ⑨ Alarm clock
- ⑩ Systolic Pressure
- ⑪ Diastolic Pressure
- ⑫ Low battery icon
- ⑬ Pulse rate
- ⑭ AF icon
- ⑮ CAM mode
- ⑯ Heartbeat icon
- ⑰ Irregular heartbeat detector
- ⑱ Memory set
- ⑲ Memory icon
- ⑳ Night time icon

- ① Day time icon
- ② Blood pressure level indicator

## Inserting the Batteries

**To insert batteries:**

1. Remove the cover of the battery compartment, which is located in the back side of the device (see the illustration).
2. Insert the batteries (4 AAA, 1.5 V), following the indicated polarity. (+ -)

### ATTENTION

- After the warning  appears, device is blocked till the batteries have been replaced.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.
- Please use «AAA» Long-Life or Alkaline 1.5V Batteries. Do not use rechargeable batteries.

## Using a Mains Adaptor

It is possible to operate Automatic Blood Pressure Monitor Model M-400A with a mains adaptor (output 5 V DC/600 mA).

1. Connect the micro plug to the socket at the device right side.
2. Plug the mains adaptor (output 5 V DC / 600 mA) to a 100-240V electrical power socket.

**NOTE:**

- Please use the adaptor of Trademark **Dr. Frei®**. Any other adaptor being used, the service center is not considered responsible for the warranty servicing/repair of the monitor.
- No power is taken from the batteries while the mains adaptor is connected to the device.
- Please consult service center if you have questions relating to the mains adaptor.

## Setting the Date and Time

---

Insert batteries or hold M/SET button for 3 seconds until the «2013» flashes on the upper right corner.

- a. When the year flashes, press M/SET button until the correct year appears. Press CAM button to confirm and change to Month setting.

01/1/2013

- b. When the month flashes, press M/SET button until the correct month appears. Press CAM button to confirm and change to date setting.

02/12/2013

- c. Set the day, hour and minute in the same way as described above.


**NOTE:** If the meter is idle for 1.5 minutes during the setting mode, it will switch off automatically.

## Setting the Alarm


---

This M-400A Series blood pressure monitor allows you to set 2 alarm clocks as a reminder to take medication or measurements.

- a. Hold M/SET then hold CAM button. Release both buttons, «1» flashes on the screen, press M/SET button to toggle between alarm «1» and «2».


- b. Press CAM button, hour flashes, press M/SET button to set the desired hours and press CAM button for the next minute setting.

- c. At the time of alarm, the meter beeps automatically. You can touch any

button to silence it immediately.

**How to switch off the alarm permanently:** hold CAM button when the alarm number «1» or «2» is flashing, release as «On» is changed to «OFF» on the display.

**NOTE:** The alarm times must be re-entered each time the batteries are replaced.

## MEASUREMENT PROCEDURE

---

### ATTENTION


- Find time to relax by sitting in a quiet atmosphere for some time before measurement.
- Efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any muscles in the arm during measurement.
- Always measure on the same arm (normally left).
- Remove any garment that fits closely to your upper arm. Do not roll the sleeve since it can squeeze your hand and this can lead to false results.
- Use only clinically approved original cuff.
- If you want to follow the results of your blood pressure measurements, always perform measurements at the same time of day, since blood pressure changes during the course of the day.
- Measurements should be done after a 5 minute rest to ensure accuracy.

### Fitting the Cuff


---

- a. Plug the cuff connector to the cuff socket on the device left side.


- b. Assemble the cuff as shown below. The smooth surface should be inside of the cuff loop and the metal D-ring should not touch your skin.

- c. Stretch your left (right) arm in front of you with your palm facing up. Slide and place the cuff onto your arm to let the air tube and artery mark region toward the lower arm. Wrap and tighten the cuff above your elbow. The artery mark on the edge of the cuff should be approx. 0.8" to 1.2" (2 cm to 3 cm) above your elbow. Align the tube over the main arteries on the inside.

- d. Leave a little free space between the arm and the cuff: you should be able to fit two fingers between them. Remove all clothing covering or constricting the measurement arm.
- e. Press the hook material firmly against the pile material. The top and bottom edges of the cuff should be tightened evenly around your upper arm.

## Taking a single measurement

---

### ATTENTION

1. Sit down for at least 10 minutes before measuring.
2. Place your elbow on a flat surface. Relax your hand with the palm facing up.
3. Make sure the cuff is about the same height as the location of your heart.

4. Remain still and do not talk or move during the measurement.


**WARNING:**

If the cuff is relatively lower (higher) than the heart, the obtained blood pressure value could be higher (lower) than the actual value.


**Always apply the pressure cuff before turning on the monitor.**

- Press START button, pump begins to inflate the cuff. The rising pressure in the cuff is shown on the display.
- Pressure in cuff increases to working pressure. Then pump stops and pressure falls.

**NOTE:** In case the inflation pressure is not sufficient to determine measurement result, the monitor automatically re-inflates to a higher level.

- The heart symbol «♥» starts to flash every heart beat once pulse is detected.
- After the measurement, systolic pressure, diastolic pressure, pulse rate and color bars classified from ESH (European Society of

- Hypertension) standards shows on the display.
- e. If irregular heartbeat was detected during the measurement, LCD displays the «» icon to remind users of heartbeat irregularity.

## **CAM (Clinically Accurate Measurement) Mode**

**Always apply the pressure cuff before turning on the monitor.**

In CAM mode, 2 measurements (at 20 sec break) are automatically taken in succession and the result is then automatically analyzed and displayed. Because people's blood pressure constantly fluctuates, a result determined in this way is more reliable than the one produced by a single measurement.


- a. Press CAM button, pump begins to inflate the cuff. LCD display shows CAM icon and the rising pressure in the cuff.


- b. Pressure in cuff increases to working pressure. Then pump stops and pressure falls.

**NOTE:** In case the inflation pressure is not sufficient to determine measurement result, the monitor automatically re-inflates to a higher level.

- c. The heart symbol «» starts to flash every heart beat once pulse is detected.


- d. After the 1<sup>st</sup> measurement, the monitor automatically starts the 2<sup>nd</sup> measurement after 20 sec countdown.

- e. After the 2<sup>nd</sup> measurement, the device analyzes the 1<sup>st</sup> and 2<sup>nd</sup> measurement values: if sys BP1- sys BP 2  $\geq$  15 mmHg, dia BP1 - dia BP2  $\geq$  10 mmHg, CAM mode will automatically start the 3<sup>rd</sup> measurement after another 20 sec countdown.

**NOTE:** In case the 1<sup>st</sup> and/or 2<sup>nd</sup> pulse rate  $\leq$  60, device will automatically take 4 times measurement in succession.

- f. The CAM measured systolic pressure, diastolic pressure, pulse rate and color bars classified from ESH (European Society of Hypertension) standards show on the display.
- g. If irregular heartbeat was detected during the measurement, LCD displays the «» icon to remind users of heartbeat irregularity.
- h. If LCD displays the «» icon, please consult your doctor to check what kinds of arrhythmia you may have. **It is important for early diagnostic of atrial fibrillation.**

### DISCONTINUING A MEASUREMENT

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the START button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

## MEMORY FUNCTION

---

### Average of All Stored Memories

---

Press M/SET button, LCD shows a set of values averaged from all the stored memories of the selected user, accompanied by an icon «A».

### Average of All Stored Day-Time Memories

---

Press M/SET button again, LCD shows another set of values averaged from all the stored day-time (4:00~9:59) taken measurements, accompanied by an icon «A» and «☀».

**NOTE:** Heart attack and stroke occurs more frequently in the morning. This feature averages the morning readings between 4:00 and 9:59AM to identify morning hypertensions.

## Average of All Stored Night-Time Memories

Press M/SET button once more, LCD shows the other set of values averaged from all the stored night-time (17:00~1:59) taken measurements, accompanied by an icon «A» and «☾».


## Viewing the Stored Values

The device stores 90 most recent blood pressure test results along with respective dates and times in the monitor memory.

### Reviewing Test Results

- Press M/SET button again after you see the average of all stored night-time memories, the single memory reading appears.
- Per the setting, the first reading you see is the last blood pressure result. Press M/SET button again, the next measurement eg. M89, as well as the previous 88 measurements (M88, M87... M1) can be displayed one after the other.


**NOTE:** CAM icon will show on LCD if the measurement is taken under CAM mode.

## Delete Memories

---

Hold down M/SET button for 7 seconds until «CL» flashes on LCD, press M/SET button again, all the stored memories inside is removed.

Before you delete all the readings stored in the memory, pls make sure you will not need to use them any more at a later time.

## ERROR MESSAGES/ TROUBLESHOOTING

---

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local distributor or **Dr. Frei**<sup>®</sup> service centre in your country. Do not attempt to repair by yourself and never try to disassemble the monitor under any circumstances.

MESSAGE	CAUSE	WHAT TO DO
	Inflation or pressure error	Refit cuff tightly and correctly.  If error still remains, contact local customer service for help.

	Blood pressure measurement error	Refit cuff tightly and correctly. Relax and repeat the measurement. If error still remains, contact local customer service for help.
	Improper movement detected during measurement (Artifact)	Relax and repeat the measurement.
	No pulse detected	Refit cuff tightly and correctly. Relax and repeat the measurement. If error still remains, contact local customer service for help.
	Problems with the monitor	Review the instructions and repeat the test.  If the monitor still does not work, please contact the local customer service for help.

## TROUBLESHOOTING

1. If no display appears after pressing START button:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.

2. If the heart rate is higher/lower than user's average:

POSSIBLE CAUSE	WHAT TO DO
Movement during measurement.	Repeat measurement.
Measurement taken just after exercise.	Rest at least 30 minutes beforehand.

3. If the result is higher/lower than user's average measurement:

POSSIBLE CAUSE	WHAT TO DO
May not be in correct position while measuring.	Adjust to the correct position to measure.
Blood pressure naturally varies from time to time.	Keep in mind for next measurement.

4. If the cuff inflates again while measuring:

POSSIBLE CAUSE	WHAT TO DO
Cuff is not fastened.	Fasten the cuff again.
If user's blood pressure is higher than the pressure the device has inflated, the device will automatically increase the pressure and start to inflate again. Stay relaxed and wait for the measurement.	

## CARE AND MAINTENANCE

---

To avoid the monitor attracting dirt, dust or other contaminants, wash and dry your hands thoroughly before use.

### Cleaning

- To clean the monitor exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do NOT flush with water.
- Do NOT use organic solvents to clean the monitor.
- Do NOT wash the pressure cuff.
- Do NOT iron the pressure cuff.

### Monitor Storage

- Storage condition: -20°C to 55°C (-4°F to 131°F), below 95% relative humidity.

- Always store or transport the monitor in its original storage case.
- Avoid dropping or heavy impact.
- Avoid direct sunlight and high humidity.

## SYMBOL INFORMATION

SYMBOL	REFERENT
	Consult instructions for use
	Manufacturer
	Serial number
	Caution, consult accompanying documents
	Dispose of in accordance with the requirements of your country
	Type BF Equipment
	Authorized representative in the European Community
	Class II
	Humidity limitation
<b>IP21</b>	Water-proof grade
	CE mark

## TECHNICAL SPECIFICATIONS

Power source	four 1.5V AAA alkaline batteries 5V 0.6A mains adaptor
Size of monitor w/o cuff	135 (L) x 95 (W) x 31.5 (H) mm
Weight	280 g with batteries
Cuff Size	M (medium): 22-38 cm (8.7-15 inches)
Memory	maximum 90 memory records
Power saving	automatic power off if system idle 1.5 min
System operating condition	10°C to 40°C (50°F to 104°F), 15-95% relative maximum humidity
Meter storage condition	-20°C to 55°C (-4°F to 131°F), 15-95% relative maximum humidity
Pressure range	20-280 mmHg
Heart rate range	40-200 beat per minute
Systolic Measurement Range	50 mmHg-250 mmHg
Diastolic Measurement Range	30 mmHg-180 mmHg
Pulse Rate Measurement Range	40-200 beats / minute
Maximum inflation pressure	280 mmHg
Accuracy of Pressure	±3 mmHg or ±2% of the reading
Accuracy of Pulse rate	±4% of the reading
Measurement unit	mmHg

This device has been tested to meet the electrical and safety requirements of:

IEC/EN 60601-1,  
IEC/EN 60601-1-2

#### Reference to Standards:

- EN 1060-1/-3, NIBP-requirements
- IEC60601-1 General requirement for safety
- IEC60601-1-2 Requirements for EMC
- EN1060-4, NIBP clinical investigation
- AAMI/ANSI/IEC 80601-2-30, ANSI/AAMI/ISO 81060-2, NIBP requirements

## WARRANTY

---

Your Automatic Blood Pressure Monitor Model M-400A is warranted for 10 years from date of purchase. Warranty for the cuff is 1 year from the date of purchase. The warranty is only valid upon presentation of the warranty card which was correctly filled in and sealed by the seller (or by the service center representative) with the date of the purchase (repair), the check, instruction manual and safe packaging of the device.

- The warranty does not apply to wear out of cuff, batteries, rubber tubes and packaging box.
- The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device.
- The device cannot be returned if it is without original packaging, with mechanical damages (scratches, stains and so on) or without any component.
- The device which was used and needs repair cannot be replaced with a new one.

In case there was violation of operating rules during warranty period, repair will be done at the expense of the consumer.

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■ Biocare Asia Co.,LTD
EC REP Shanghai International Trading Corp. GmbH (Hamburg)
Manufactured by the order of: Medhouse Swiss GmbH,